**OsteoBiol<sup>®</sup>**  
by Tecnos

## Bone Grafting Materials

REGENERATION SCIENCE

INSPIRED BY NATURE

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- (1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H  
**PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE**  
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19
- (2) NANNMARK U, SENNERBY L  
**THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS**  
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70. EPUB 2008 JAN 30
- (3) CRESPI R, CAPPARÈ P, ROMANOS GE, MARIANI E, BENASCIUTTI E, GHERLONE E  
**CORTICOCANCELLOUS PORCINE BONE IN THE HEALING OF HUMAN EXTRACTION SOCKETS: COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST GENE EXPRESSION PROFILES IN VIVO**  
INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG; 26(4):866-72
- (4) CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G  
**BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS**  
CLIN ORAL IMPLANTS RES, 2014 MAY 26 EPUB AHEAD OF PRINT
- (5) FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M  
**COMPARISON OF A XENOGENIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE**  
MATER SCI ENG C MATER BIOL APPL, 2013 AUG 1;33(6):3506-13
- (6) FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S  
**INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING - A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION**  
CLIN ORAL IMPLANTS RES, 2014 SEP 15 EPUB AHEAD OF PRINT
- (7) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R  
**PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS**  
CLIN IMPLANT DENT RELAT RES, 2011 NOV 14, EPUB AHEAD OF PRINT
- (8) CARDAROPOLI D, CARDAROPOLI G  
**PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY**  
INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77
- (9) BARONE A, CRESPI R, NICOLI ALDINI N, FINI M, GIARDINO R, COVANI U  
**MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS**  
INT J ORAL MAXILLOFAC IMPLANTS, 2005, 20: 519-525
- (10) BOTTINI LP, RICCI L, PIATTELLI A, PERROTTI V, IEZZI G  
**BUCCO-LINGUAL CRESTAL BONE CHANGES AROUND IMPLANTS IMMEDIATELY PLACED IN FRESH SOCKETS IN ASSOCIATION OR NOT WITH PORCINE BONE: A NON-BLINDED RANDOMIZED CONTROLLED TRIAL IN HUMANS**  
J PERIODONTOL, 2012 OCT 29, EPUB AHEAD OF PRINT

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- (1) BARONE A, RICCI M, COVANI U, NANNMARK U, AZARMEHR I, CALVO GUIRADO JL  
MAXILLARY SINUS AUGMENTATION USING PREHYDRATED CORTICOCANCELLOUS PORCINE BONE: HISTOMORPHOMETRIC EVALUATION AFTER 6 MONTHS  
CLIN IMPLANT DENT RELAT RES, 2012 JUN;14(3):373-9 Epub 2010 MAY 11
- (2) NANNMARK U, SENNERBY L  
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS  
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70. Epub 2008 JAN 30
- (3) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C  
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAPTS USED IN MAXILLARY SINUS FLOOR ELEVATION  
CLIN ORAL IMPLANTS RES, 2012 JAN 26, Epub AHEAD OF PRINT
- (4) CALVO GUIRADO JL, RAMIREZ FERNANDEZ MP, NEGRI B, DELGADO RUIZ RA, MATÉ SANCHEZ DE VAL JE, GOMEZ MORENO G  
EXPERIMENTAL MODEL OF BONE RESPONSE TO COLLAGENIZED XENOGRAPTS OF PORCINE ORIGIN (OSTEOBIOL® mp3®): A RADIOLOGICAL AND HISTOMORPHOMETRIC STUDY  
CLIN IMPLANT DENT RELAT RES, 2013 FEB;15(1):143-51 Epub 2011 MAR 31
- (5) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, DELGADO RUIZ RA, MATÉ SANCHEZ DE VAL JE, VICENTE ORTEGA V, MESEGUER OLMOS L  
BONE RESPONSE TO HYDROXYAPATITES WITH OPEN POROSITY OF ANIMAL ORIGIN (PORCINE OSTEOBIOL® mp3®) AND BOVINE (ENDOBON®): A RADIOLOGICAL AND HISTOMORPHOMETRIC STUDY  
CLIN ORAL IMPLANTS RES, 2011 JUL; 22(7):767-73
- (6) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H  
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE  
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19
- (7) SILVESTRI M, MARTEGANI P, D'AVENIA F, FARNETI M, CAPRI D, PAOLANTONI G, LANDI L  
SIMULTANEOUS SINUS AUGMENTATION WITH IMPLANT PLACEMENT: HISTOMORPHOMETRIC COMPARISON OF TWO DIFFERENT GRAFTING MATERIALS. A MULTICENTER DOUBLE-BLIND PROSPECTIVE RANDOMIZED CONTROLLED CLINICAL TRIAL  
INT J ORAL MAXILLOFAC IMPLANTS, 2013 MAR-APR;28(2):543-9
- (8) PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A, VERROCCHI D, VOLPE S, SENNERBY L  
A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOSS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY  
CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58 Epub 2010 OCT 26
- (9) BARONE A, ALDINI NN, FINI M, GIARDINO R, CALVO GUIRADO JL, COVANI U  
XENOGRAFT VERSUS EXTRACTION ALONE FOR RIDGE PRESERVATION AFTER TOOTH REMOVAL: A CLINICAL AND HISTOMORPHOMETRIC STUDY  
J PERIODONTOL, 2008 AUG;79(8):1370-7
- (10) BARONE A, COVANI U  
MAXILLARY ALVEOLAR RIDGE RECONSTRUCTION WITH NONVASCULARIZED AUTOGENOUS BLOCK BONE: CLINICAL RESULTS  
J ORAL MAXILLOFAC SURG, 2007, 65:2039-2046
- (11) WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T  
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES  
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7

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- (1) CALVO GUIRADO JL, GÓMEZ MORENO G, GUARDIA J, ORTIZ RUIZ A, PIATTELLI A, BARONE A, MARTÍNEZ GONZÁLEZ JM, MESENGUER OLMO L, LÓPEZ MARÍ L, DORADO CB  
**BIOLOGICAL RESPONSE TO PORCINE XENOGRIFT IMPLANTS: AN EXPERIMENTAL STUDY IN RABBITS**  
IMPLANT DENT, 2012 APR;21(2):112-7
- (2) ARCURI C, CECCHETTI F, GERMANO F, MOTTA A, SANTACROCE C  
**CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS**  
MINERVA STOMATOL, 2005, 54: 351-362.
- (3) NANNMARK U, AZARMEHR I  
**SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECT**  
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3
- (4) BARONE A, AMERI S, COVANI U  
**IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS**  
EUR J IMPLANT PROSTHODONTICS, 2006, 2: 99-106
- (5) CASSETTA M, RICCI L, IEZZI G, DELLAQUILA D, PIATTELLI A, PERROTTI V  
**RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN**  
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (6) CALVO GUIRADO JL, PARDO ZAMORA G, SAEZ YUGUERO MR  
**RIDGE SPLITTING TECHNIQUE IN ATROPHIC ANTERIOR MAXILLA WITH IMMEDIATE IMPLANTS, BONE REGENERATION AND IMMEDIATE TEMPORISATION: A CASE REPORT**  
JIR DENT ASSOC, 2007 WINTER, 53(4):187-90

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- (1) BARONE A, CORNELINI R, CIAGLIA R, COVANI U  
IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND  
SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE  
SERIES  
INT J PERIODONTICS RESTORATIVE DENT, 2008 JUN; 28(3):283-9
- (2) SANTAGATA M, GUARINIELLO L, RAUSO R, TARTARO G  
IMMEDIATE LOADING OF DENTAL IMPLANT AFTER SINUS FLOOR  
ELEVATION WITH OSTEOTOME TECHNIQUE: A CLINICAL REPORT  
AND PRELIMINARY RADIOGRAPHIC RESULTS  
J ORAL IMPLANTOL, 2010 DEC; 36(6):485-489
- (3) COVANI U, CORNELINI R, BARONE A  
BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS  
WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH  
INT J ORAL MAXILLOFAC IMPLANTS, 2008 SEP-OCT; 23(5):841-6
- (4) CARDAROPOLI D, CARDAROPOLI G  
HEALING OF GINGIVAL RECESIONS USING A COLLAGEN  
MEMBRANE WITH A DEMINERALIZED XENOGRFT: A RANDOMIZED  
CONTROLLED CLINICAL TRIAL  
INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB;29(1):59-67
- (5) NANNMARK U, AZARMEHR I  
SHORT COMMUNICATION: COLLAGENATED CORTICOCANCEL-  
LOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY  
DEFECTS  
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3

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- (1) TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S  
**THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS**  
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91
- (2) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R  
**OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS**  
EUR J INFLAMM, 2011, VOL. 9, NO. 3(S), 103-107
- (3) ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A  
**HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL**  
J PERIODONTOL 2006 DEC;77(12):1984-90
- (4) SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G  
**MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS**  
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18
- (5) EZZI G, DEGIDI M, PIATTELLI A, MANGANO C, SCARANO A, SHIBLI JA, PERROTTI V  
**COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS**  
CLIN ORAL IMPLANTS RES, 2012 DEC;23(12):1369-76 Epub 2011 NOV 2
- (6) BARONE A, AMERI S, COVANI U  
**IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS**  
EUR J IMPLANT PROSTHODONTICS, 2006;2: 99-106
- (7) BARONE A, TOTI P, QUARANTA A, DERCHI G, COVANI U  
**THE CLINICAL OUTCOMES OF IMMEDIATE VERSUS DELAYED RESTORATION PROCEDURES ON IMMEDIATE IMPLANTS: A COMPARATIVE COHORT STUDY FOR SINGLE-TOOTH REPLACEMENT**  
CLIN IMPLANT DENT RELAT RES, 2014 MAY 22 Epub AHEAD OF PRINT
- (8) BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, CALVO GUIRADO JL, NEGRI B, DI FELICE R, COVANI U  
**VOLUMETRIC ANALYSIS OF REMODELLING PATTERN AFTER RIDGE PRESERVATION COMPARING USE OF TWO TYPES OF XENOGRAFTS. A MULTICENTRE RANDOMIZED CLINICAL TRIAL**  
CLIN ORAL IMPLANTS RES, 2015 FEB 26 Epub AHEAD OF PRINT
- (9) SCARANO A, PIATTELLI A, ASSENZA B, QUARANTA A, PERROTTI V, PIATTELLI M, IEZZI G  
**PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION**  
J ORAL MAXILLOFAC SURG, 2010 AUG; 68(8):1869-73

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(2) FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M  
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT  
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66

(3) ESPOSITO M, CANNIZZARO G, SOARDI E, PISTILLI R, PIATTELLI M, CORVINO V, FELICE P  
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM-LONG, 4 MM-WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL  
EUR J ORAL IMPLANTOL, 2012;5(1):19-33

(4) FELICE P, PISTILLI R, PIATTELLI M, SOARDI E, CORVINO V, ESPOSITO M  
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A RANDOMISED CONTROLLED TRIAL  
EUR J ORAL IMPLANTOL, SUMMER, 5(2):149-61

(5) PISTILLI R, FELICE P, PIATTELLI M, GESSAROLI M, SOARDI E, BARAUSSE C, BUTI J, CORVINO V, ESPOSITO M  
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. ONE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL  
EUR J ORAL IMPLANTOL, 2013;6(4):343-357

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- (1) NANNMARK U, SENNERBY L  
**THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS**  
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70
- (2) CARDAROPOLI D, CARDAROPOLI G  
**HEALING OF GINGIVAL RECESIONS USING A COLLAGEN MEMBRANE WITH A DEMINERALIZED XENOGRFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL**  
INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB; 29(1):59-67
- (3) BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G  
**FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL**  
CLIN ORAL IMPLANTS RES, 2014 MAR 1 EPUB AHEAD OF PRINT
- (4) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U  
**TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING**  
CLIN ORAL IMPLANTS RES, 2012 JUL 12, EPUB AHEAD OF PRINT
- (5) CARDAROPOLI D, CARDAROPOLI G  
**PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY**  
INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77
- (6) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C  
**ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAFTS USED IN MAXILLARY SINUS FLOOR ELEVATION**  
CLIN ORAL IMPLANTS RES, 2012 JAN 26, EPUB AHEAD OF PRINT
- (7) BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U  
**A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL**  
CLIN ORAL IMPLANTS RES, 2013 JAN;24(1):1-6. EPUB 2011 DEC 12
- (8) SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G  
**MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS**  
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18
- (9) CASSETTA M, RICCI L, IEZZI G, CALASSO S, PIATTELLI A, PERROTTI V  
**USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES**  
INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC;32(6):E182-8
- (10) COVANI U, MARCONCINI S, CRESPI R, BARONE A  
**IMMEDIATE IMPLANT PLACEMENT AFTER REMOVAL OF A FAILED IMPLANT: A CLINICAL AND HISTOLOGICAL CASE REPORT**  
J ORAL IMPLANTOL, 2009; 35(4):189-95
- (11) FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M  
**VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT**  
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66

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- (1) RINNA C, REALE G, FORESTA E, MUSTAZZA MC  
**ORBITAL WALL RECONSTRUCTION WITH SWINE BONE CORTEX**  
J CRANIOFAC SURG, 2009 MAY; 20(3): 881-4
- (2) PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A, VERROCCHI D, VOLPE S, SENNERBY L  
**A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOSS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY**  
CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58. EPUB 2010 OCT 26
- (3) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R  
**PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS**  
CLIN IMPLANT DENT RELAT RES, 2011 NOV 14, EPUB AHEAD OF PRINT
- (4) WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T  
**THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES**  
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7
- (5) SCARANO A, PIATTELLI M, CARINCI F, PERROTTI V  
**REMOVAL, AFTER 7 YEARS, OF AN IMPLANT DISPLACED INTO THE MAXILLARY SINUS. A CLINICAL AND HISTOLOGIC CASE REPORT**  
J OF OSSEOINTEGRATION, 2009
- (6) HINZE M, VRIELINCK L, THALMAIR T, WACHTEL H, BOLZ W  
**ZYGOMATIC IMPLANT PLACEMENT IN CONJUNCTION WITH SINUS BONE GRAFTING: THE "EXTENDED SINUS ELEVATION TECHNIQUE". A CASE-COHORT STUDY**  
ORAL CRANIOFAC TISSUE ENG, 2011;1:188-197
- (7) SCARANO A, PIATTELLI A, IEZZI G, VARVARA G  
**SPONTANEOUS BONE FORMATION ON THE MAXILLARY SINUS FLOOR IN ASSOCIATION WITH SURGERY TO REMOVE A MIGRATED DENTAL IMPLANT: A CASE REPORT**  
MINERVA STOMATOL, 2014 OCT;63(10):351-9
- (8) RINNA C, UNGARI C, SALTAREL A, CASSONI A, REALE G  
**ORBITAL FLOOR RESTORATION**  
J CRANIOFAC SURG, 2005, 16(6): 968-972
- (9) OZEL B, FINDIKCIOGLU K, SEZGIN B, GUNAY K, BARUT I, OZMEN S  
**A NEW OPTION FOR THE RECONSTRUCTION OF ORBITAL FLOOR DEFECTS WITH HETEROLOGOUS CORTICAL BONE**  
J CRANIOMAXILLOFAC SURG, 2015 OCT;43(8):1583-8 EPUB 2015 JUL 4

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[www.bbb.education](http://www.bbb.education)

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# INNOVACIÓN

A close-up photograph of a scientist in a white lab coat, hairnet, and face mask, looking through a microscope. The background is a blurred laboratory setting with various glassware and equipment. The word "INNOVACIÓN" is overlaid in large, bold, blue capital letters on the left side of the image.

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- (1) ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A  
**HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL**  
J PERIODONTOL, 2006 DEC;77(12):1984-90
- (2) SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G  
**MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS**  
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18
- (3) IEZZI G, DEGIDI M, PIATTELLI A, MANGANO C, SCARANO A, SHIBLI JA, PERROTTI V  
**COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS**  
CLIN ORAL IMPLANTS RES, 2012 DEC;23(12):1369-76. EPUB 2011 NOV 2
- (4) CALVO GUIRADO JL, GOMEZ MORENO G, GUARDIA J, ORTIZ RUIZ A, PIATTELLI A, BARONE A, MARTÍNEZ GONZÁLEZ JM, MESENGUER OLMO L, LÓPEZ MARÍ L, DORADO CB  
**BIOLOGICAL RESPONSE TO PORCINE XENOGRFT IMPLANTS: AN EXPERIMENTAL STUDY IN RABBITS**  
IMPLANT DENT, 2012 APR;21(2):112-7
- (5) FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M  
**COMPARISON OF A XENOGENEIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE**  
MATER SCI ENG C MATER BIOL APPL, 2013 AUG 1;33(6):3506-13. EPUB 2013 MAY 3
- (6) BARONE A, CRESPI R, ALDINI NN, FINI M, GIARDINO R, COVANI U  
**MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS**  
INT J ORAL MAXILLOFAC IMPLANTS, 2005 JUL-AUG; 20(4):519-25
- (7) CRESPI R, CAPPARÈ P, ROMANOS GE, MARIANI E, BENASCIUTTI E, GHERLONE E  
**CORTICOCANCELLOUS PORCINE BONE IN THE HEALING OF HUMAN EXTRACTION SOCKETS: COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST GENE EXPRESSION PROFILES IN VIVO**  
INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG; 26(4):866-72
- (8) CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G  
**BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS**  
CLIN ORAL IMPLANTS RES, 2014 MAY 26 EPUB AHEAD OF PRINT

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(1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H  
**PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE**  
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19

(2) TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S  
**THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS**  
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91

(3) NANNMARK U, SENNERBY L  
**THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS**  
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70. EPUB 2008 JAN 30

(4) CRESPI R, CAPPARÈ P, ROMANOS GE, MARIANI E, BENASCIUTTI E, GHERLONE E  
**CORTICOCANCELLOUS PORCINE BONE IN THE HEALING OF HUMAN EXTRACTION SOCKETS: COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST GENE EXPRESSION PROFILES IN VIVO**  
INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG; 26(4):866-72

(5) FBARONE A, RICCI M, COVANI U, NANNMARK U, AZARMEHR I, CALVO GUIRADO JL  
**MAXILLARY SINUS AUGMENTATION USING PREHYDRATED CORTICOCANCELLOUS PORCINE BONE: HISTOMORPHOMETRIC EVALUATION AFTER 6 MONTHS**  
CLIN IMPLANT DENT RELAT RES, 2012 JUN;14(3):373-9. EPUB 2010 MAY 11

(6) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R  
**OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS**  
EUR J INFLAMMAT, 2011, VOL. 9, NO. 3(S), 103-107

(7) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U  
**TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING**  
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7. EPUB 2012 JUL 12

(8) BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G  
**FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL**  
CLIN ORAL IMPLANTS RES, 2014 MAR 1 EPUB AHEAD OF PRINT

(9) FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S  
**INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING - A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION**  
CLIN ORAL IMPLANTS RES, 2014 SEP 15 EPUB AHEAD OF PRINT

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(1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H  
**PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE**  
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19

(2) BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U  
**A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL**  
CLIN ORAL IMPLANTS RES, 2013 JAN;24(1):1-6. EPUB 2011 DEC 12

(3) SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G  
**MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS**  
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18

(4) CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G  
**BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS**  
CLIN ORAL IMPLANTS RES, 2014 MAY 26 EPUB AHEAD OF PRINT

(5) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R  
**OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS**  
EUR J INFLAMMAT, 2011, VOL. 9, NO. 3(S), 103-107

(6) HSU FY, CHUEH SC, WANG YJ  
**MICROSPHERES OF HYDROXYAPATITE/RECONSTITUTED COLLAGEN AS SUPPORTS FOR OSTEOBLAST CELL GROWTH**  
BIOMATERIALS 1999, 20:1931-1936

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
(1) GRIFFITH LG, NAUGHTON G  
**TISSUE ENGINEERING-CURRENT CHALLENGES AND EXPANDING OPPORTUNITIES**  
 SCIENCE 2002, 295:1009-14

(2) REDDI AH  
**MORPHOGENESIS AND TISSUE ENGINEERING OF BONE AND CARTILAGE: INDUCTIVE SIGNALS, STEM CELLS, AND BIOMIMETIC BIOMATERIALS**  
 TISSUE ENG 2000, 6(4):351-59

(3) NAKASHIMA N, REDDI AH  
**THE APPLICATION OF BONE MORPHOGENETIC PROTEINS TO DENTAL TISSUE ENGINEERING**  
 NAT BIOTECHNOL 2003, 9:1025-32

(4) SALASZNYK RM, WILLIAMS WA, BOSKEY A, BATORSKY A, PLOPPER GE  
**ADHESION TO VITRONECTIN AND COLLAGEN I PROMOTES OSTEOGENIC DIFFERENTIATION OF HUMAN MESENCHYMAL STEM CELLS**  
 J BIOMED BIOTECHNOL 2004, 1:24-34

(5) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R  
**OSTEOBIOL<sup>®</sup> INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS**  
 EUR J INFLAMMAT, 2011, VOL. 9, NO. 3(S), 103-107

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(1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H  
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE  
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19

(2) NANNMARK U, SENNERBY L  
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS  
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70. EPUB 2008 JAN 30

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**iDisfrute de [osteobiol.com](http://osteobiol.com)!**

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- 1 | COVANI U, AMERI S, CRESPI R, BARONE A  
**PRESERVAZIONE DEL PROCESSO ALVEOLARE CON OSSO ESTEREOLOGO. CONSIDERAZIONI ISTOLOGICHE**  
ITALIAN ORAL SURGERY, 2004, VOL 3, 1: 17-23
- 2 | CASSETTA M, CALASSO S, VOZZA I, DELL'AQUILA D  
**REHABILITATION OF ATROPHIC ALVEOLAR CRESTS WITH CYLINDRICAL SANDBLASTED AND ACID ETCHED IMPLANTS: A PILOT STUDY**  
EUR J IMPLANT PROSTHODONTICS, 2005;(3)1:133-144
- 3 | ARCURI C, CECCHETTI F, GERMANO F, MOTTA A, SANTACROCE C  
**CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC BONE SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS**  
MINERVA STOMATOL, 2005 JUN;54(6):351-62
- 4 | BARONE A, CRESPI R, ALDINI NN, FINI M, GIARDINO R, COVANI U  
**MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS**  
INT J ORAL MAXILLOFAC IMPLANTS, 2005 JUL-AUG; 20(4):519-25
- 5 | RINNA C, UNGARI C, SALTARELLA, CASSONI A, REALE G  
**ORBITAL FLOOR RESTORATION**  
J CRANIOFAC SURG, 2005 NOV; 16(6):968-72
- 6 | BARONE A, AMERI S, COVANI U  
**IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS**  
EUR J IMPLANT PROSTHODONTICS, 2006,2: 99-106
- 7 | BARONE A, SANTINI S, SBORDONE L, CRESPI R, COVANI U  
**A CLINICAL STUDY OF THE OUTCOMES AND COMPLICATIONS ASSOCIATED WITH MAXILLARY SINUS AUGMENTATION**  
INT J ORAL MAXILLOFAC IMPLANTS, 2006 JAN-FEB; 21(1):81-5
- 8 | COVANI U, BARONE A, CORNELINI R, CRESPI R  
**CLINICAL OUTCOME OF IMPLANTS PLACED IMMEDIATELY AFTER IMPLANT REMOVAL**  
J PERIODONTOL, 2006 APR;77(4):722-7
- 9 | ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A  
**HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL**  
J PERIODONTOL, 2006 DEC;77(12):1984-90
- 10 | TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S  
**THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS**  
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91
- 11 | BARONE A, COVANI U  
**MAXILLARY ALVEOLAR RIDGE RECONSTRUCTION WITH NON- VASCULARIZED AUTOGENOUS BLOCK BONE: CLINICAL RESULTS**  
J ORAL MAXILLOFAC SURG, 2007 OCT;65(10):2039-46
- 12 | CALVO GUIRADO JL, PARDO ZAMORA G, SAEZ YUGUERO MR  
**RIDGE SPLITTING TECHNIQUE IN ATROPHIC ANTERIOR MAXILLA WITH IMMEDIATE IMPLANTS, BONE REGENERATION AND IMMEDIATE TEMPORISATION: A CASE REPORT**  
J IR DENT ASSOC 2007 WINTER;53(4):187-90
- 13 | DEL CORSO M  
**SOFT TISSUE RESPONSE TO PLATELET RICH FIBRIN: CLINICAL EVIDENCES**  
COSMETIC DENTISTRY, 2008, 3:16-20
- 14 | BARONE A, SANTINI S, MARCONCINI S, GIACOMELLI L, GHERLONE E, COVANI U  
**OSTEOTOMY AND MEMBRANE ELEVATION DURING THE MAXILLARY SINUS AUGMENTATION PROCEDURE. A COMPARATIVE STUDY: PIEZOELECTRIC DEVICE VS. CONVENTIONAL ROTATIVE INSTRUMENTS**  
CLIN ORAL IMPLANTS RES, 2008 MAY;19(5):511-5.
- 15 | BARONE A, CORNELINI R, CIAGLIA R, COVANI U  
**IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE SERIES**  
INT J PERIODONTICS RESTORATIVE DENT, 2008 JUN; 28(3):283-9
- 16 | BARONE A, ALDINI NN, FINI M, GIARDINO R, CALVO GUIRADO JL, COVANI U  
**XENOGRAFT VERSUS EXTRACTION ALONE FOR RIDGE PRESERVATION AFTER TOOTH REMOVAL: A CLINICAL AND HISTOMORPHOMETRIC STUDY**  
J PERIODONTOL, 2008 AUG;79(8):1370-7
- 17 | COVANI U, CORNELINI R, BARONE A  
**BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH**  
INT J ORAL MAXILLOFAC IMPLANTS, 2008 SEP-OCT; 23(5):841-6
- 18 | CARDAROPOLI D, CARDAROPOLI G  
**PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY**  
INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77
- 19 | NANNMARK U, SENNERBY L  
**THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS**  
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70.
- 22 | CALVO GUIRADO JL, GOMEZ MORENO G, BARONE A, CUTANDO A, ALCARAZ BANOS M, CHIVA F, LOPEZ MARI L, GUARDIA J  
**MELATONIN PLUS PORCINE BONE ON DISCRETE CALCIUM DEPOSIT IMPLANT SURFACE STIMULATES OSTEOINTEGRATION IN DENTAL IMPLANTS**  
J PINEAL RES, 2009, 47(2):164-72
- 23 | SCARANO A, PIATTELLI M, CARINCI F, PERROTTI V  
**REMOVAL, AFTER 7 YEARS, OF AN IMPLANT DISPLACED INTO THE MAXILLARY SINUS. A CLINICAL AND HISTOLOGIC CASE REPORT**  
J OSSEOINTEGR, 2009;1(1):35-40
- 24 | COVANI U, MARCONCINI S, CRESPI R, BARONE A  
**IMMEDIATE IMPLANT PLACEMENT AFTER REMOVAL OF A FAILED IMPLANT: A CLINICAL AND HISTOLOGICAL CASE REPORT**  
J ORAL IMPLANTOL, 2009; 35(4):189-95
- 25 | CALVO GUIRADO JL, GOMEZ MORENO G, LOPEZ MARI L, ORTIZ RUIZ AJ, GUARDIA J  
**ATRAUMATIC MAXILLARY SINUS ELEVATION USING THREADED BONE DILATORS FOR IMMEDIATE IMPLANTS. A THREE-YEAR CLINICAL STUDY**  
MEDICINA ORAL, PATOLOGIA ORAL Y CIRUGIA BUCAL, 2010 MAR 1; 15(2):E366-70
- 26 | FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H  
**PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE**  
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19
- 27 | GRENGA PL, REALE G, COFONE C, MEDURI A, CERUTI P, GRENGA R  
**HESS AREA RATIO AND DIPLOPIA: EVALUATION OF 30 PATIENTS UNDERGOING SURGICAL REPAIR FOR ORBITAL BLOW-OUT FRACTURE**  
OPHTHAL PLAST RECONSTR SURG, 2009 MAR-APR; 25(2):123-5
- 28 | CRESPI R, CAPPARÈ P, GHERLONE E  
**DENTAL IMPLANTS PLACED IN EXTRACTION SITES GRAFTED WITH DIFFERENT BONE SUBSTITUTES: RADIOGRAPHIC EVALUATION AT 24 MONTHS**  
J PERIODONTOL, 2009 OCT; 80(10):1616-1621
- 29 | RINNA C, REALE G, FORESTA E, MUSTAZZA MC  
**MEDIAL ORBITAL WALL RECONSTRUCTION WITH SWINE BONE CORTEX**  
J CRANIOFAC SURG, 2009 MAY; 20(3): 881-4
- 30 | CARDAROPOLI D, CARDAROPOLI G  
**HEALING OF GINGIVAL RECESSIONS USING A COLLAGEN MEMBRANE WITH A THE MINERALIZED XENOGRAFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL**  
INT J OF PERIODONTICS AND RESTORATIVE DENTISTRY, 2009 FEB; 29(1):59-67
- 31 | NANNMARK U, AZARMEHR I  
**SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS**  
CLINICAL IMPLANT DENTISTRY AND RELATED RESEARCH, 2010 JUN 1; 12(2):161-3
- 32 | BARONE A, RICCI M, COVANI U, NANNMARK U, AZARMEHR I, CALVO GUIRADO JL  
**MAXILLARY SINUS AUGMENTATION USING PRE-HYDRATED CORTICOCANCELLOUS PORCINE BONE: HISTOMORPHOMETRIC EVALUATION AFTER 6 MONTHS**  
CLIN IMPLANT DENT RELAT RES, 2012 JUN;14(3):373-9
- 33 | CALVO GUIRADO JL, GOMEZ MORENO G, LOPEZ MARI L, GUARDIA J, MARINEZ GONZALEZ JM, TRESGUERRRES IF, PAREDES SD, FUENTES BRETO L  
**ACTIONS OF MELATONIN MIXED WITH COLLAGENIZED PORCINE BONE VERSUS PORCINE BONE ONLY ON OSTEOINTEGRATION OF DENTAL IMPLANTS**  
J PINEAL RES, 2010, 48:194-203
- 34 | SCARANO A, PIATTELLI A, ASSENZA B, QUARANTA A, PERROTTI V, PIATTELLI M, IEZZI G  
**PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION**  
J ORAL MAXILLOFAC SURG, 2010 AUG; 68(8):1869-73
- 35 | ROSSI R, MORALES RS, FRASCARIA M, BENZI R, SQUADRITO N  
**PLANNING IMPLANTS IN THE ESTHETIC ZONE USING A NEW IMPLANT 3D NAVIGATION SYSTEM**  
EUR J ESTHETIC DENT, 2010 SUMMER; 5(2):172-88
- 37 | SCARANO A, CARINCI F, ASSENZA B, PIATTELLI M, MURMURA G, PIATTELLI A  
**VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE USING AN INLAY TECHNIQUE WITH A XENOGRAFT WITHOUT MINISCREWS AND MINIPLATES: CASE SERIES**  
CLIN ORAL IMPLANTS RES, 2011 OCT;22(10):1125-30
- 38 | PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A, VERROCCHI D, VOLPE S, SENNERBY L  
**A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY**  
CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58
- 39 | SANTAGATA M, GUARINIELLO L, TARTARO G  
**A MODIFIED EDENTULOUS RIDGE EXPANSION (MERE) TECHNIQUE FOR IMMEDIATE PLACEMENT OF IMPLANTS. A CASE REPORT**  
J ORAL IMPLANTOL, 2011 MAR;37 SPEC N°:114-9
- 40 | BARONE A, RICCI M, CALVO GUIRADO JL, COVANI U  
**BONE REMODELLING AFTER REGENERATIVE PROCEDURES AROUND IMPLANTS PLACED IN FRESH EXTRACTION SOCKETS: AN EXPERIMENTAL STUDY IN BEAGLE DOGS**  
CLIN ORAL IMPLANTS RES, 2011 OCT; 22(10):1131-7
- 41 | SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G  
**MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRIC EVALUATION AFTER 4 AND 6 MONTHS**  
CLIN IMPLANT DENT RELAT RES, 2011 MAR;13(1):13-18
- 42 | CALVO GUIRADO JL, RAMIREZ FERNANDEZ MP, NEGRI B, DELGADO RUIZ RA, MATÉ SANCHEZ DE VAL JE, GOMEZ MORENO G  
**EXPERIMENTAL MODEL OF BONE RESPONSE TO COLLAGENIZED XENOGRAFTS OF PORCINE ORIGIN (OSTEOBIOL® MP3®): A RADIOLOGICAL AND HISTOMORPHOMETRIC STUDY**  
CLIN IMPLANT DENT RELAT RES, 2013 FEB;15(1):143-51
- 43 | RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, DELGADO RUIZ RA, MATÉ SANCHEZ DE VAL JE, VICENTE ORTEGA V, MESEGUER OLMOS L  
**BONE RESPONSE TO HYDROXYAPATITES WITH OPEN POROSITY OF ANIMAL ORIGIN (PORCINE OSTEOBIOL® MP3®) AND BOVINE (ENDOBON®): A RADIOLOGICAL AND HISTOMORPHOMETRIC STUDY**  
CLIN ORAL IMPLANTS RES, 2011 JUL; 22(7):767-73

- 44 | CRESPI R, CAPPARÈ P, ROMANOS GE, MARIANI E, BENASCIUTTI E, GHERLONE E  
**CORTICOCANCELOUS PORCINE BONE IN THE HEALING OF HUMAN EXTRACTION SOCKETS: COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST GENE EXPRESSION PROFILES IN VIVO**  
INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG; 26(4):866-72
- 45 | HINZE M, VRIELINCK L, THALMAIR T, WACHTEL H, BOLZ W  
**ZYGOMATIC IMPLANT PLACEMENT IN CONJUNCTION WITH SINUS BONE GRAFTING: THE "EXTENDED SINUS ELEVATION TECHNIQUE". A CASE-COHORT STUDY**  
ORAL CRANIOFAC TISSUE ENG, 2011; 1:188-197
- 46 | IEZZI G, DEGIDI M, PIATTELLI A, MANGANO C, SCARANO A, SHIBLI JA, PERROTTI V  
**COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS**  
CLIN ORAL IMPLANTS RES, 2012 DEC;23(12):1369-76
- 47 | SLOTTÉ C, LINDFORS N, NANNMARK U  
**SURGICAL RECONSTRUCTION OF PERI-IMPLANT BONE DEFECTS WITH PREHYDRATED AND COLLAGENATED PORCINE BONE AND COLLAGEN BARRIERS: CASE PRESENTATIONS**  
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):714-23
- 48 | BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U  
**A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL**  
CLIN ORAL IMPLANTS RES, 2013 JAN; 24(1):1-6
- 49 | SANTAGATA M, GUARINIELLO L, RAUSO R, TARTARO G  
**IMMEDIATE LOADING OF DENTAL IMPLANT AFTER SINUS FLOOR ELEVATION WITH OSTEOTOME TECHNIQUE: A CLINICAL REPORT AND PRELIMINARY RADIOGRAPHIC RESULTS**  
J ORAL IMPLANTOL, 2010 DEC; 36(6):485-489
- 50 | FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R  
**PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS**  
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):707-13
- 51 | RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C  
**ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAFTS USED IN MAXILLARY SINUS FLOOR ELEVATION**  
CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30
- 52 | CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PERROTTI V  
**RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN**  
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- 53 | BARONE A, ORLANDO B, CINGANO L, MARCONCINI S, DERCHI G, COVANI U  
**A RANDOMIZED CLINICAL TRIAL TO EVALUATE AND COMPARE IMPLANTS PLACED IN AUGMENTED VS. NON-AUGMENTED EXTRACTION SOCKETS. A 3-YEAR EVALUATION**  
J PERIODONTOL, 2012 JUL;83(7):836-46
- 54 | CALVO GUIRADO JL, MATÉ SANCHEZ JE, DELGADO RUIZ L, RAMIREZ FERNANDEZ MP  
**CALCULATION OF BONE GRAFT VOLUME USING 3D RECONSTRUCTION SYSTEM**  
MED ORAL PATOL CIR BUCAL, 2011 MAR 1, 16(2):E260-4
- 55 | ESPOSITO M, CANNIZZARO G, SOARDI E, PISTILLI R, PIATTELLI M, CORVINO V, FELICE P  
**POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM-LONG, 4 MM-WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL**  
EUR J ORAL IMPLANTOL, 2012 SPRING;5(1):19-33
- 56 | FELICE P, PIANA L, CHECCHI L, PISTILLI R, PELLEGRINO G  
**VERTICAL RIDGE AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH A 2-STAGE INLAY TECHNIQUE: A CASE REPORT**  
IMPLANT DENT, 2012 JUN;21(3):190-5
- 57 | BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U  
**TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING**  
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7
- 58 | CASSETTA M, RICCI L, IEZZI G, CALASSO S, PIATTELLI A, PERROTTI V  
**USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES**  
INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC;32(6):E182-8
- 59 | BOTTINI LP, RICCI L, PIATTELLI A, PERROTTI V, IEZZI G  
**BUCCO-LINGUAL CRESTAL BONE CHANGES AROUND IMPLANTS IMMEDIATELY PLACED IN FRESH SOCKETS IN ASSOCIATION OR NOT WITH PORCINE BONE: A NON-BLINDED RANDOMIZED CONTROLLED TRIAL IN HUMANS**  
J PERIODONTOL, 2012 OCT 29, EPUB AHEAD OF PRINT
- 60 | BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R  
**OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS**  
EUR J INFLAMMAT, 2011, VOL. 9, NO. 3 (5), 103-107
- 61 | FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M  
**VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT**  
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66
- 62 | FICKL S, JOCKEL-SCHNEIDER Y, LINCKE T, BECHTOLD M, FISCHER KR, SCHLAGENHAUF U  
**PORCINE DERMAL MATRIX FOR COVERING OF RECESSION TYPE DEFECTS: A CASE SERIES**  
QUINTESSENCE INT, 2013;44(3):243-6
- 63 | SILVESTRI M, MARTEGANI P, D'AVENIA F, FARNETI M, CAPRI D, PAOLANTONI G, LANDI L  
**SIMULTANEOUS SINUS AUGMENTATION WITH IMPLANT PLACEMENT: HISTOMORPHOMETRIC COMPARISON OF TWO DIFFERENT GRAFTING MATERIALS. A MULTICENTER DOUBLE-BLIND PROSPECTIVE RANDOMIZED CONTROLLED CLINICAL TRIAL**  
INT J ORAL MAXILLOFAC IMPLANTS, 2013 MAR-APR;28(2):543-9
- 64 | WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T  
**THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES**  
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7
- 65 | RODRIGUEZ JG, ELDBANY RM  
**VERTICAL SPLITTING OF THE MANDIBULAR BODY AS AN ALTERNATIVE TO INFERIOR ALVEOLAR NERVE LATERALIZATION**  
INT J ORAL MAXILLOFAC SURG, 2013 SEP;42(9):1060-6
- 66 | FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M  
**COMPARISON OF A XENOGENEIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE**  
MATER SCI ENG C, MATER BIOL APP, 2013 AUG 1;33(6):3506-13
- 67 | FELICE P, PISTILLI R, PIATTELLI M, SOARDI E, CORVINO V, ESPOSITO M  
**POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A RANDOMISED CONTROLLED TRIAL**  
EUR J ORAL IMPLANTOL, SUMMER;5(2):149-61
- 68 | TRAINI T, PIATTELLI A, CAPUTI S, DEGIDI M, MANGANO C, SCARANO A, PERROTTI V, IEZZI G  
**REGENERATION OF HUMAN BONE USING DIFFERENT BONE SUBSTITUTE BIOMATERIALS**  
CLIN IMPLANT DENT RELAT RES, 2013 MAY 17
- 69 | CALVO GUIRADO JL, GOMEZ MORENO G, GUARDIA J, ORTIZ RUIZ A, PIATTELLI A, BARONE A, MARTÍNEZ GONZÁLEZ JM, MESENGUER OLMO L, LÓPEZ MARÍ L, DORADO CB  
**BIOLOGICAL RESPONSE TO PORCINE XENOGRAFT IMPLANTS: AN EXPERIMENTAL STUDY IN RABBITS**  
IMPLANT DENT, 2012 APR;21(2):112-7
- 70 | KOLMAS J, SZWAJA M, KOŁODZIEJSKI W  
**SOLID-STATE NMR AND IR CHARACTERIZATION OF COMMERCIAL XENOGENEIC BIOMATERIALS USED AS BONE SUBSTITUTES**  
J PHARM BIOMED ANAL, 2012 MAR 5;61:136-41
- 71 | PISTILLI R, FELICE P, PIATTELLI M, GESSAROLI M, SOARDI E, BARAUSSE C, BUTI J, CORVINO V, ESPOSITO M  
**POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. ONE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL**  
EUR J ORAL IMPLANTOL, 2013 WINTER;6(4):343-357
- 72 | FICKL S, NANNMARK U, SCHLAGENHAUF U, HÜRZELER M, KEBSCHULL M  
**PORCINE DERMAL MATRIX IN THE TREATMENT OF DEHISCENCE-TYPE DEFECTS – AN EXPERIMENTAL SPLIT-MOUTH ANIMAL TRIAL**  
CLIN ORAL IMPLANTS RES, 2014 FEB 19
- 73 | BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G  
**FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL**  
CLIN ORAL IMPLANTS RES, 2014 MAR 1
- 74 | BARONE A, RICCI M, ROMANOS GE, TONELLI P, ALFONSI F, COVANI U  
**BUCCAL BONE DEFICIENCY IN FRESH EXTRACTION SOCKETS: A PROSPECTIVE SINGLE COHORT STUDY**  
CLIN ORAL IMPLANTS RES, 2014 MAR 31
- 75 | BARONE A, TOTI P, QUARANTA A, DERCHI G, COVANI U  
**THE CLINICAL OUTCOMES OF IMMEDIATE VERSUS DELAYED RESTORATION PROCEDURES ON IMMEDIATE IMPLANTS: A COMPARATIVE COHORT STUDY FOR SINGLE-TOOTH REPLACEMENT**  
CLIN IMPLANT DENT RELAT RES, 2014 MAY 22
- 76 | CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G  
**BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS**  
CLIN ORAL IMPLANTS RES, 2014 MAY 26
- 77 | GHENO E, PALERMO A, BUFFOLI B, RODELLA LF  
**THE EFFECTIVENESS OF THE USE OF XENOGENEIC BONE BLOCKS MIXED WITH AUTOLOGOUS CONCENTRATED GROWTH FACTORS (CGF) IN BONE REGENERATION TECHNIQUES**  
J OSSEOINTEGRATION 2014;6(2):37-42
- 78 | FALISI G, GALLI M, VITTORINI-VELASQUEZ P, GALLEGOS-RIVERA JC, MINASI R, DE BIASI A, DI PAOLO C  
**USE OF 3D CARTILAGE SCAFFOLDS FOR THE STABILIZATION OF IMPLANTS AND BONE REGENERATION WITH THE FIT-LOCK TECHNIQUE**  
ACTA ODONTOL LATINOAM 2013;26(3):167-172
- 79 | FISCHER KR, FICKL S, MARDAS N, BOZEC L, DONOS N  
**STAGE-TWO SURGERY USING COLLAGEN SOFT TISSUE GRAFTS: CLINICAL CASES AND ULTRASTRUCTURAL ANALYSIS**  
QUINTESSENCE INT, 2014 NOV-DEC; 45(10):853-60
- 80 | FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S  
**INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING – A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION**  
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1135-42.

81 | SCARANO A, MURMURA G, SINIARI B, ASSENZA B, SOLLAZZO V, SPINELLI G, CARINCI F

**EXPANSION OF THE ALVEOLAR BONE CREST WITH ULTRASONIC SURGERY DEVICE: CLINICAL STUDY IN MANDIBLE**

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**DELAYED EXPANSION OF THE ATROPHIC MANDIBLE BY ULTRASONIC SURGERY: A CLINICAL AND HISTOLOGIC CASE SERIES**

INT J ORAL MAXILLOFAC IMPLANTS, 2015 JAN-FEB;30(1):144-9

83 | LORENZON G, BUTTARELLO GM, CHESSA G

**CASE REPORT: IMPLANT PLACEMENT AND IMMEDIATE LOADING WITH SIMULTANEOUS BONE REGENERATION FOLLOWING JAW ODONTOGENIC CYST ENUCLEATION DENTISTRY, 2015, 5:2**

84 | THALWAIR T, FICKL S, SCHNEIDER D, HINZE M, WACHTEL H  
**DIMENSIONAL ALTERATIONS OF EXTRACTION SITES AFTER DIFFERENT ALVEOLAR RIDGE PRESERVATION TECHNIQUES – A VOLUMETRIC STUDY**

J CLIN PERIODONTOL, 2013 JUL;40(7):721-7

85 | MANESCU A, GIULIANI A, MOHAMMADI S, TROMBA G, MAZZONI S, DIOMEDE F, ZINI N, PIATELLI A, TRUBIANI O  
**OSTEOGENIC POTENTIAL OF DUAL-BLOCKS CULTURED WITH HUMAN PERIODONTAL LIGAMENT STEM CELLS: IN VITRO AND SYNCHROTRON**

J PERIODONTAL RES, 2015 JUN 11

86 | SCARANO A, PIATELLI A, IEZZI G, VARVARA G  
**SPONTANEOUS BONE FORMATION ON THE MAXILLARY SINUS FLOOR IN ASSOCIATION WITH SURGERY TO REMOVE A MIGRATED DENTAL IMPLANT: A CASE REPORT**

MINERVA STOMATOL, 2014 OCT;63(10):351-9

87 | BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, CALVO GUIRADO JL, NEGRİ B, DI FELICE R, COVANI U  
**VOLUMETRIC ANALYSIS OF REMODELLING PATTERN AFTER RIDGE PRESERVATION COMPARING USE OF TWO TYPES OF XENOGRAFTS. A MULTICENTRE RANDOMIZED CLINICAL TRIAL**

CLIN ORAL IMPLANTS RES, 2015 FEB 26

88 | ESPOSITO M, GRUSOVIN MG, LAMBERT F, MATOS S, PIETRUSKA M, ROSSI R, SALHI L, BUTI J  
**THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRABONY DEFECT - A MULTICENTRE RANDOMISED CONTROLLED TRIAL**

EUR J ORAL IMPLANTOL, 2015;8(3):233-244

89 | OZEL B, FINKICIOGLU K, SEZGIN B, GUNEY K, BARUT I, OZMEN S

**A NEW OPTION FOR THE RECONSTRUCTION OF ORBITAL FLOOR DEFECTS WITH HETEROLOGOUS CORTICAL BONE**

J CRANIO-MAXILLOFAC SURG, 2015 OCT;43(8):1583-8 EPUB 2015 JUL 4

90 | CALVO GUIRADO JL, AGUILAR SALVATIERRA A, RAMIREZ FERNANDEZ M, MATÉ SANCHEZ DEL VAL JE, DELGADO RUIZ RA, GOMEZ MORENO G  
**BONE RESPONSE TO COLLAGENATED XENOGRAFTS OF PORCINE ORIGIN (MP3®) AND A BOVINE BONE MINERAL GRAFTING (4BONETM XBM) GRAFTS IN TIBIA DEFECTS:**

**EXPERIMENTAL STUDY IN RABBITS**

CLIN ORAL IMPLANTS RES, 2015 OCT;43(8):1583-8.

91 | CORBELLA S, TASCHERI S, WEINSTEIN R, DEL FABBRO M  
**HISTOMORPHOMETRIC OUTCOMES AFTER LATERAL SINUS FLOOR ELEVATION PROCEDURE: A SYSTEMATIC REVIEW OF THE LITERATURE AND META-ANALYSIS**

CLIN ORAL IMPLANTS RES, 2015 OCT 10

92 | CALVO GUIRADO JL, GOMEZ MORENO G, MATÉ SANCHEZ JE, LOPEZ MARI L, DELGADO RUIZ R, ROMANOS GE  
**NEW BONE FORMATION IN BONE DEFECTS AFTER MELATONIN AND PORCINE BONE**

**GRAFTS: EXPERIMENTAL STUDY IN RABBITS**

CLIN ORAL IMPLANTS RES, 2015 APR;26(4):399-406.

93 | BARONE A, MARCONCINI S, GIAMMARINARO E, MUIRITSKY E, GELPI F, COVANI U

**CLINICAL OUTCOMES OF IMPLANTS PLACED IN EXTRACTION SOCKETS AND IMMEDIATELY RESTORED: A 7-YEAR SINGLE-COHORT PROSPECTIVE STUDY**

CLIN IMPLANT DENT RELAT RES, 2016 FEB 16.

94 | ÇAKIR M, KARACA IR, AŞEĞÜL F, KAYMAZ F, BOZKAYA S  
**EXPERIMENTAL EVALUATION OF THE EFFECTS OF ANKAFERD BLOOD STOPPER AND COLLAGENATED HETEROLOGOUS BONE GRAFT ON BONE HEALING IN SINUS FLOOR AUGMENTATION**

CLIN ORAL IMPLANTS RES, 2015 MAR-APR;30(2):279-85

95 | LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F  
**REGENERATION OF ATROPHIC CRESTAL RIDGES WITH RESORBABLE LAMINA: TECHNICAL NOTE**

J BIOL REGUL HOMEOST AGENTS 2015 JUL-SEP;29(3 SUPPL 1):97-100

96 | ETTORRE V, DE MARCO P, ZARA S, FERROTTI V, SCARANO A, DI CRESCENZO A, PETRINI M, HADAD C, BOSCO D, ZAVAN B, VALBONETTI L, SPOTO G, IEZZI G, PIATELLI A, CATALDI A, FONTANA A

**IN VITRO AND IN VIVO CHARACTERIZATION OF GRAPHENE OXIDE COATED PORCINE BONE GRANULES**

CARBON, JULY 2016, VOLUME 103, PAGES 291–298 AVAILABLE ONLINE 14 MARCH 2016

97 | ROSSI R, RANCITELLI D, POLI PP, RASIA DAL POLO M, NANINMARK U, MAIORANA C

**THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY**

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**BONE REGENERATION IN ILIAC CRESTAL DEFECTS: AN EXPERIMENTAL STUDY ON SHEEP**

BIOMED RESEARCH INTERNATIONAL 2016

99 | FELICE P, ZUCCHELLI G, CANNIZZARO G, BARAUSSÉ C, DIAZZI M, TRULLENQUE-ERIKSSON A, ESPOSITO M

**IMMEDIATE, IMMEDIATE-DELAYED (6 WEEKS) AND DELAYED (4 MONTHS) POST-EXTRACTIVE SINGLE IMPLANTS: 4-MONTH POST-LOADING DATA FROM A RANDOMISED CONTROLLED TRIAL**

EUROPEAN JOURNAL OF ORAL IMPLANTOLOGY

2016;9(3):233-247

100 | ROMBOUTS C, JEANNEAU C, CAMILLERI J, LAURENT P,

**ABOUT CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENEIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS**

DENTAL MATERIALS JOURNAL 2016 DEC 1;35(6):900-907

101 | BARONE A, TOTI P, MARCONCINI S, DERCHI G, MARCHIONNI S, COVANI U

**ESTHETIC OUTCOME OF IMPLANTS PLACED IN FRESH EXTRACTION SOCKETS BY CLINICIANS WITH OR WITHOUT EXPERIENCE: A MEDIUM-TERM RETROSPECTIVE EVALUATION**

INT JOURNAL OF ORAL AND MAXILLOFACIAL IMPLANTS

2016 NOV/DIC;31(6):1397-1406

102 | TALLARICO M, XHANARI E, PISANO M, DE RIU G, TULLIO A, MELONI SM

**SINGLE POST-EXTRACTIVE ULTRA-WIDE 7 MM-DIAMETER IMPLANTS VERSUS IMPLANTS PLACED IN MOLAR HEALED SITES AFTER SOCKET PRESERVATION FOR MOLAR REPLACEMENT: 6-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL**

EUROPEAN JOURNAL OF ORAL IMPLANTOLOGY

2016;9(3):263-275

103 | LOPEZ MA, MANZULLI N, CASALE M, ORMANER Z, CARINCI F  
**THE USE OF RESORBABLE HETEROLOGOUS CORTICAL LAMINA AS A NEW SINUS LIFT FLOOR: A TECHNICAL NOTE**

JOURNAL OF BIOLOGICAL REGULATORS AND HOMEOSTATIC AGENTS 2016 APR-JUN;30(2 SUPPL 1):75-79

104 | LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F, ORMANER Z, LAURTANO D

**THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES**

JOURNAL OF BIOLOGICAL REGULATORS AND HOMEOSTATIC AGENTS 2016 APR-JUN;30(2 SUPPL 1):81-85

105 | ESPOSITO M, ZUCCHELLI G, BARAUSSÉ C, PISTILLI R, TRULLENQUE-ERIKSSON A, FELICE P

**FOUR MM-LONG VERSUS LONGER IMPLANTS IN AUGMENTED BONE IN POSTERIOR ATROPHIC INJAYS: 4-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL**

EUROPEAN JOURNAL OF ORAL IMPLANTOLOGY

2016;9(4):393-409

106 | BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, NEGRİ B, DI FELICE R, MARCHIONNI S, CALVO-GUIRADO JL, COVANI U, NANINMARK U

**CLINICAL AND HISTOLOGICAL CHANGES AFTER RIDGE PRESERVATION WITH TWO XENOGRAFTS: PRELIMINARY RESULTS FROM A MULTICENTER RANDOMIZED CONTROLLED CLINICAL TRIAL**

JOURNAL OF CLINICAL PERIODONTOLOGY 2017

FEB;44(2):204-214

107 | ALFONSI F, BORGIA V, IEZZI G, PIATELLI A, COVANI U, TONELLI P, BARONE A

**MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING THE MANAGEMENT OF EXTRACTION SOCKETS**

CURRENT PHARMACEUTICAL BIOTECHNOLOGY 2017;18(1):64-75

108 | BARONE A, TOTI P, MENCHINI FABRIS GB, MARCHIONNI S, COVANI U

**EARLY VOLUMETRIC CHANGES AFTER VERTICAL**

**AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH INTERPOSITIONAL BLOCK GRAFT VERSUS ONLAY BONE GRAFT: A RETROSPECTIVE RADIOLOGICAL STUDY**

JOURNAL OF CRANIO-MAXILLOFACIAL SURGERY 2017 JAN 25

109 | BARONE A, TOTI P, FUNEL N, CAMPANI D, COVANI U  
**EXPRESSION OF SP7, RUNX1, DLX5, AND CTNINB1 IN HUMAN MESENCHYMAL STEM CELLS CULTURED ON XENOGENEIC BONE SUBSTITUTE AS COMPARED WITH MACHINED TITANIUM**

IMPLANT DENTISTRY 2014 AUG;23(4):407-15

110 | LOZANO CARRASCAL N, DELGADO RUIZ RA, GARGALLO ALBIOL J, MATÉ SANCHEZ JE, HERNANDEZ ALFARO F, CALVO GUIRADO JL

**XENOGRAFTS SUPPLEMENTED WITH PAMINDRONATE PLACED IN POSTEXTRACTION SOCKETS TO AVOID CRESTAL BONE RESORPTION. EXPERIMENTAL STUDY IN FOX HOUND DOGS**

CLINICAL ORAL IMPLANTS RESEARCH 2016 FEB;27(2):149-55

111 | SLOTTE C, NANINMARK U

**AUGMENTATION OF SEVERLY RESORBED ALVEOLAR CRESTS WITH BOVINE BONE BLOCKS: A HUMAN HISTOLOGY PILOT STUDY**

CLINICAL IMPLANT DENTISTRY AND RELATED RESEARCH 2017 MAR 1

112 | ESPOSITO M, ZUCCHELLI G, CANNIZZARO G, CHECCHI L, BARAUSSÉ C, TRULLENQUE-ERIKSSON A, FELICE P

**IMMEDIATE, IMMEDIATE-DELAYED (6 WEEKS) AND DELAYED (4 MONTHS) POST-EXTRACTIVE SINGLE IMPLANTS: 1-YEAR POST-LOADING DATA FROM A RANDOMISED CONTROLLED TRIAL**

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**BONE REGENERATION INDUCED BY BONE PORCINE BLOCK WITH BONE MARROW STROMAL STEM CELLS IN A MINIPIG MODEL OF MANDIBULAR “CRITICAL SIZE” DEFECT**

STEM CELLS INTERNATIONAL 2017;2017:9082869

114 | SCARANO A

**TRADITIONAL POSTEXTRACTIVE IMPLANT SITE PREPARATION COMPARED WITH PRE-EXTRACTIVE INTERRADICULAR IMPLANT BED PREPARATION IN THE MANDIBULAR MOLAR REGION, USING AN ULTRASONIC DEVICE: A RANDOMIZED PILOT STUDY**

INT JOURNAL OF ORAL AND MAXILLOFACIAL IMPLANTS

2017 MAY/JUN;32(3):655-660

115 | BARONE A, TOTI P, MENCHINI FABRIS GB, DERCHI G, MARCONCINI S, COVANI U

**EXTRA ORAL DIGITAL SCANNING AND IMAGING SUPERIMPOSITION FOR VOLUME ANALYSIS OF BONE REMODELING AFTER TOOTH EXTRACTION WITH AND WITHOUT 2 TYPES OF PARTICULATE PORCINE MINERAL INSERTION: A RANDOMIZED CONTROLLED TRIAL**

CLINICAL IMPLANT DENTISTRY AND RELATED RESEARCH 2017 MAY 16

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